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EN. SURGICAL LIGHT HANDLE COVER

INDICATIONS

Surgical light handle covers are a medical device intended to prevent medical staff intentionally or accidentally touching the handle of the light during a surgical procedure. The special toothed configuration of the inner plastic disk makes this device adaptable to any handle with diameter between 20 mm and 40 mm. Different lengths are available compatible with different handle measurements.

REF	MEASUREMENTS
LHC-01	Light Handle Cover (12 cm pouch length)
LHC-03	Light Handle Cover (15 cm pouch length)

INSTRUCTIONS FOR USE

After opening the package, remove the product from the inner pouch.
Insert the end of the handle into the protective cover by pressing the plastic support.
Move the plastic support up to the top of the handle so that the handle is completely covered.

DURATION OF USE

For single use only. The product can be used until the end of the surgical intervention.

WARNINGS

This product must be used by authorised medical staff.
Once in place, do not reposition. If the product is repositioned, the fastening system may work loose and the sheath may become detached in the surgical field.
Follow the instructions for use and disposal of the product.
Product for single use only, sterilized by ethylene oxide. Do not re-use. Do not re-sterilize.
The re-use or re-sterilization of this product could cause:
- Cross infection due to the presence of biological waste;
- Alterations to materials;
- Loss of functionality of the product.
Re-sterilization of the product is not guaranteed.
Do not use if the individual packaging is open or damaged.
After having checked the integrity of the packaging, verify its content. If the packaging shows signs of damage or visible defects, do not use the product and return it to TELIC, S.A.U.
Check the expiry date on the packaging. Do not use after the expiry date.
The recommended storage and usage conditions are between 5° C and 40° C (41° F - 104° F). The relative humidity restrictions are 20-80%.

DISPOSAL

- The product must be disposed of according to hospital protocols for contaminated waste, in accordance with the applicable regulations.
- Dispose of the packaging in accordance with recyclable waste management.

WARRANTY

TELIC, S.A.U. guarantees that the product complies with the Regulation (EU) 2017/745 on Medical devices and was manufactured following the procedures of the TELIC, S.A.U. quality management system certified under ISO standard 13485. TELIC, S.A.U. accepts no liability for any medical costs or direct or indirect harm due to the lack of operation or incorrect operation of the products when used without observing the instructions for use. We advise you to duly report any operational problem or product fault to the Quality Assurance Department of TELIC, S.A.U. In the event of any serious incident related to the product, this must be reported immediately to TELIC, S.A.U. (incidents@telic.es) and/or the competent authority of the Member State in which the incident occurs.

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Emergo Consulting (UK) Limited