TELIC, S.A.U. guarantees that the product complies with the Regulation (EU) 2017/745 on Medical devices and was manufactured following the procedures of the TELIC, S.A.U. quality management system certified under ISO standard 13485. TELIC, S.A.U. accepts no liability for any medical costs or direct or indirect harm due to the lack of operation or incorrect operation of the products when used without observing the instructions or the instructions for use. We advise you to duly report any operational problem or product fault to the Quality Assurance Department of TELIC, S.A.U. in the event of any serious incident related to the product, in which the incident occurs.

**WARRANTY**  
TELIC, S.A.U. guarantees that the product complies with the Regulation (EU) 2017/745 on Medical devices and was manufactured following the procedures of the TELIC, S.A.U. quality management system certified under ISO standard 13485. TELIC, S.A.U. accepts no liability for any medical costs or direct or indirect harm due to the lack of operation or incorrect operation of the products when used without observing the instructions or the instructions for use. We advise you to duly report any operational problem or product fault to the Quality Assurance Department of TELIC, S.A.U. in the event of any serious incident related to the product, in which the incident occurs.  
- This unit must be reported immediately to TELIC, S.A.U. (incidents@telics.es) and/or the competent authority of the Member State  
- Recycle packaging  
- Once you have removed the plastic sheath from the electrodes, dispose of it immediately in the appropriate waste container.  
- Dispose off as hospital contaminated waste, in accordance with regulations in force.  
**DISPOSAL**  
- Once you have removed the plastic sheath from the electrodes, dispose of it immediately in the appropriate waste container.  
- Recycle packaging  
- Dispose off as hospital contaminated waste, in accordance with regulations in force.  
- During surgical procedures electrodes should be placed as far as possible from the elector- surgical area to minimize unwarranted RF current flow and burnings in the patient.  
- Electrodes may not be repositioned.  
- In the event of use over an extended period, check for secure adhesion and absence of skin reaction on a daily basis  
- Use electrodes, never allow contact with an external electrical source.  
- Do not use the electrodes if the pouch is found damaged.

**PRECAUTIONS**  
- Do not use the electrodes if the pouch is found damaged.  
- Do not use the gel if it is dry or discoloured.  
- Dissolve the adhesive.  
- Careless removal of electrode may cause damage to the skin. Should the electrode be difficult to remove, use water to remove the adhesive.  
- Gentle removal of electrodes from the pouch to be ensured so as to prevent any injury to the fingers as plastic covers are not compatible with conventional cleaning and sterilization procedures.  
- Procedures may leave harmful residues or cause malfunction of the device. The design and materials used are not suitable for reuse.  
- The product is not for internal use. In case of accidental ingestion of the electrode consult a physician immediately.  
- For single use only. Use on other patients can cause cross infection. Do not soak, reuse, or sterilize this device, as these procedures may leave harmful residues or cause malfunction of the device.  
- The product is not for internal use. In case of accidental ingestion of the electrode consult a physician immediately.  
- And use only on intact, clean skin (do not use on open wounds, lesions, infected or inflamed areas).  
- Any allergic or other skin reactions should be immediately evaluated by a physician.  
- ECG electrodes should only be used by or in consultation with a health care provider familiar with their proper placement

**WARNINGS**  
- ECG electrodes should only be used by or in consultation with a health care provider familiar with their proper placement  
- If they are stored in this way, the electrodes can be used safely up to 7 days after the bag has been opened.  
- Once the bag has been opened, use the electrodes as soon as possible.  
- Keep the remaining electrodes in a closed plastic container, protected from moisture and within the indicated temperature margins.  
**STORAGE**  
ECG Electrodes may be used on patients for up to 24 hours.

**DURATION OF USE**  
- Remove gel residue, if required, with a paper towel or water.  
- After monitoring, lift the tab to gently remove the electrode from the skin. (5)  
- Turn de equipment on and perform the ECG procedure. (4)  
- Place the electrodes on the skin. Smooth the backing surface with the fingers as necessary. (3)  
- Connect the lead wires to the electrodes and peel the plastic liner. (2)  
- application of the electrodes. (1)  
- If necessary, clean the skin with appropriate preparation material. Always ensure the application site is dry prior to application.  
- Determine the application sites. Remove any excessive hair, calluses and ensure that the skin is clear of any oils, creams or fluids, dirt or other impediments to good electrode contact.  
- Open the pouch along the indicated line.  
**AREA FOR APPLICATION / CONTRAINDICATIONS**  
- Connect the lead wires to the electrodes and peel the plastic liner. (2)  
- application of the electrodes. (1)  
- If necessary, clean the skin with appropriate preparation material. Always ensure the application site is dry prior to application.  
- Determine the application sites. Remove any excessive hair, calluses and ensure that the skin is clear of any oils, creams or fluids, dirt or other impediments to good electrode contact.  
- Open the pouch along the indicated line.

**INSTRUCTIONS FOR USE**  
- Due to high adhesion strength of the adhesive material, this product is not suitable for use on patients with fragile skin. Not suitable for use on newborn or premature infants.  
- For use in passive sensing of bio-electric signals from the intact skin.  
**INTENDED USE**  
ECG Electrodes may be used on patients for up to 24 hours.

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**PRECAUTIONS**  
- Do not use the gel if it is dry or discoloured.  
- Dissolve the adhesive.  
- Careless removal of electrode may cause damage to the skin. Should the electrode be difficult to remove, use water to remove the adhesive.  
- Gentle removal of electrodes from the pouch to be ensured so as to prevent any injury to the fingers as plastic covers are not compatible with conventional cleaning and sterilization procedures.  
- Procedures may leave harmful residues or cause malfunction of the device. The design and materials used are not suitable for reuse.  
- The product is not for internal use. In case of accidental ingestion of the electrode consult a physician immediately.  
- For single use only. Use on other patients can cause cross infection. Do not soak, reuse, or sterilize this device, as these procedures may leave harmful residues or cause malfunction of the device.  
- The product is not for internal use. In case of accidental ingestion of the electrode consult a physician immediately.  
- And use only on intact, clean skin (do not use on open wounds, lesions, infected or inflamed areas).  
- Any allergic or other skin reactions should be immediately evaluated by a physician.  
- ECG electrodes should only be used by or in consultation with a health care provider familiar with their proper placement

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**INSTRUCTIONS FOR USE**  
- Due to high adhesion strength of the adhesive material, this product is not suitable for use on patients with fragile skin. Not suitable for use on newborn or premature infants.  
- For use in passive sensing of bio-electric signals from the intact skin.  
**INTENDED USE**  
ECG Electrodes may be used on patients for up to 24 hours.  
**DURATION OF USE**  
- Remove gel residue, if required, with a paper towel or water.  
- After monitoring, lift the tab to gently remove the electrode from the skin. (5)  
- Turn de equipment on and perform the ECG procedure. (4)  
- Place the electrodes on the skin. Smooth the backing surface with the fingers as necessary. (3)  
- Connect the lead wires to the electrodes and peel the plastic liner. (2)  
- application of the electrodes. (1)  
- If necessary, clean the skin with appropriate preparation material. Always ensure the application site is dry prior to application.  
- Determine the application sites. Remove any excessive hair, calluses and ensure that the skin is clear of any oils, creams or fluids, dirt or other impediments to good electrode contact.  
- Open the pouch along the indicated line.

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