EN. SET OF TWO DEFIBRILLATION ELECTRODES FOR INFANT/CHILD

INSTRUCTIONS FOR USE

1. Clean and dry the application site. Avoid contact with alcohol or other disinfectants.
2. Open the package and remove the electrodes.
3. Remove the protective film from electrodes.
4. Place the electrodes on the patient at clinically indicated site, making sure that all adhesive area is stuck to the skin. Be careful not to trap any pockets of air between the hydrogel and the patient's skin.
5. Connect the electrodes to the defibrillator or multi-function cable.
6. Follow the defibrillator instructions regarding the required function.

USE TIME

In case that repeated defibrillations are required, electrodes may be left on the patient during several hours depending on the skin conditions. Replace the electrode pads after 24 hours on skin.

WARNINGS

Verify connector type. — Do not open package until immediately prior to use — Do not use if package is open or damaged — Do not use if electrode pad is damaged or hydrogel is dry — Operator should be thoroughly familiar with instructions of use prior to use — This product

may be used by authorized medical personnel only or for accredited personnel for use of Automatic External Defibrillators (AEDs) — Do not

fold, bend or puncture the defibrillation electrodes — Follow the defibrillator instructions — Do not use defibrillator settings higher than 4 J/kg — Keep the defibrillation electrodes away from other electrodes or metallic parts in contact with the patient — Verify that no one is

in contact with the patient or anything that touches the patient while defibrillating. (This includes defibrillation electrodes, cables, bed, etc.) — Electromedical equipment not being protected against defibrillation shock should be disconnected from patient before defibrillating — Misuse or misapplication of defibrillation electrodes may result in patient burns or ineffective therapy — Effective defibrillation may cause visible reddening under the application area of defibrillation electrodes on the patient skin — Recommended limits of use between 5°C (41°F) and 35°C (95°F) — Non sterile product. Do not sterilize — Possible explosion and fire hazard if used in presence of flammable anaesthetics or concentrated oxygen — Check the expiry date on the defibrillation electrodes package — Avoid defibrillation electrodes placement near the generator of an internal pacemaker — Place electrodes according to operating instructions of defibrillator or pacemaker manufacturer — Figures are guidelines for positioning defibrillation electrodes and are not intended to replace clinical preferences. — If possible do not apply electrode pads on damaged skin — In case of patient movement it is necessary to check defibrillation electrodes adhesion — Defibrillation electrodes are not repositionable. Replace with new electrode pads if repositioning is required — Apply a new pair of defibrillation electrodes when there is an indication of poor pad adhesion — Remove electrode pads by slowly peeling from patient's skin — Discard after use, as per hospital's routine — Reuse of this product may cause failures in the defibrillation and/or cross-contamination to the patient.

DISPOSAL

The product must be disposed according to hospital protocols for contaminated waste, in accordance with regulations in force.

— Plastic liner should be immediately disposed off in the assigned waste bin for disposal after removal from the electrodes.

— Recycle package.

WARRANTY

TELIC, S.A.U. guarantees that the product complies with Directive 93/42/EEC on Medical devices and was manufactured following the procedures

of the TELIC, S.A.U. quality management system certified under ISO

standard 13485. TELIC, S.A.U. accepts no liability for any medical costs

or direct or indirect harm due to the lack of operation or incorrect

operation of the products when used without observing the instructions

for use. We advise you to duly report any operational problem or product

fault to the Quality Assurance Department of TELIC, S.A.U. in the event

of any serious incident related to the product, this must be reported

immediately to TELIC, S.A.U. (incidents@telic.es) and/or the competent

authority of the Member State in which the incident occurs.

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