# blayco®

## ES. MANGO DE BISTURÍ ELECTROQUIRÚRGICO

### INFORMACIÓN IMPORTANTE

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## EN. ELECTROSURGICAL PENCIL

### IMPORTANT INFORMATION

Please carefully read these instructions before use.

### DEVICE DESCRIPTION

Electrosurgical pencils are used for the cutting and coagulation of soft tissue in electrosurgery.

The product is provided with a conductive cable which is designed for connection to the high frequency electrosurgical generator (Electrosurgical unit). The device is intended to be used in combination with a neutral electrode to close the circuit.

BLAYCO® Electrosurgical pencils can be used in combination with any of BLAYCO® Electrosurgical accessories.

The maximum applicable voltage of this product is 5kVp.

### INTENDED USE

Electrosurgical pencil is intended to use for cutting and coagulation to remove tissue and control bleeding, by using high frequency current supplied by the electrosurgical generator, during electrosurgical procedures. The medical device is disposable for single use and supplied as sterile with an electrode tip.

### INDICATION

Electrosurgical pencils are indicated for use in tissue cutting and coagulating as required or encountered in General surgery, thoracic surgery, gynecology surgery provided that the operating area is easily accessible by the tip of the electrosurgical pencil.

### MODELS SPECIFICATION

This medical device includes several references as per table below. Please refer to the primary label for the identification of product

Reference	Description	Tip Shape	Exposed length	Cable length	Connector/ plug
MB-100	Hand control electrosurgical pencil	70 mm Blade	4 cm	3 m	standard 3 pin
MB-100/5	Hand control electrosurgical pencil	70 mm Blade	4 cm	5 m	standard 3 pin
MB-300	Hand control electrosurgical pencil	70mm insulated blade	4 cm	3 m	standard 3 pin
MB-400	Rocker switch electrosurgical pencil	70mm insulated blade	4 cm	3 m	standard 3 pin
MB-500	Rocker switch electrosurgical pencil	70 mm Blade	4 cm	3 m	standard 3 pin

### CONTRAINDICATIONS

The product is contraindicated in cases where the patients have a cardiac pacemaker or other active implants, there is a potential risk of interferences or damages.

Electrosurgical pencils are contraindicated for use in patients with automatic implantable cardioverter/ defibrillator (AICD).

Also, the use of this device is contraindicated in patients with the following conditions: Acute inflammation, pregnancy, malignancy, unresolved/ unaccessible adnexal pathological. This product is contraindicated for use in tubal ligation procedure.

### PRECAUTIONS

1. Examine pencil and electrode for defects prior to use.
2. Apparent low output signal or failure of the high-frequency electrosurgical equipment to function correctly at the normal operating settings, may indicate a faulty application of the neutral electrode or poor contact in connections. In this case, the application of the neutral electrode and connections should be checked before selecting a higher output power.
3. When more than one active electrodes are plugged into same generator, be aware activation of one pencil may take the receptacle of all the other pencils "live". Read and understand the warnings and cautions given by the electrosurgical generator manufacturer.
4. Do not prepare the operation site with flammable solutions, such as alcohol.
5. This device is for single patient use only. It is not intended for disinfection and subsequent reuse, which may result in device failure and/or risk of contamination.
6. This device is not intended for reprocessing or re-sterilization. Reprocessing and/or re-sterilization may damage the product, rendering it unusable and/or may lead to device failure, which could result inpatient illness, injury or death.
7. Visibly exposed metal of the active electrode shaft where it connects with the active handle.

### INSTRUCTIONS FOR USE

1. The electrosurgical pencil is designed for use with electrosurgical generators equipped with the standard three-pin connection. Prior to use, carefully read the manufacturer instructions on how to operate the electrosurgical unit and accessories.
2. Visually examine the package for damage before opening the package. Replace the pencil if the package or the pencil itself shows signs of damage.
3. Check the production and expiry date printed on the package. In case that the product is expired, replace by a backup device within expire date.
4. Open the package using aseptic technique.
5. Uncoil the cable. Examine the pencil and cable for damage. Replace the pencil if any sign of damage is found.

- Remove and discard the electrode tip cap. Ensure that the electrode is securely seated in the pencil.
- Plug the cable connector firmly into the corresponding plug of the electrosurgical generator.
- Verify that the "cut" control on the pencil works when the CUT button (yellow) or switch is pressed. Repeat the same for the coagulation control COAG (blue button or switch).
- Discard the pencil after use as a biohazard. Do not re-sterilize or reuse.

#### WARNINGS

- Electrosurgery is dangerous. Misuses of any part of the electrosurgical system may result in serious patient burns or injuries. Read and understand all instructions for use and warnings before using the product.
- Operator must ensure that connected accessories can be used with at least the maximum peak output voltage of the high-frequency electrosurgical equipment that can be set at the intended output control setting with the intended operating mode, with reference to the diagrams required.
- The entire area of the neutral electrode should be reliably adhered to the patient's body and placed as close to the operating field as possible.
- The patient should not come into contact with metal parts which are earthed or which have an appreciable capacitance to earth (for example operating table supports, etc.). The use of antistatic sheeting is recommended for this purpose.
- Skin-to-skin contact (for example between the arms and body of the patient) should be avoided, for example by the insertion of dry gauze.
- When high-frequency electrosurgical equipment and physiological monitoring equipment are simultaneously used on the same patient, any monitoring electrodes should be placed as far as possible from the electrosurgical electrodes. Needle monitoring electrodes are not recommended.
- The cables of the electrosurgical electrodes should be positioned in such a way that contact with the patient or other leads is avoided. Temporarily unused active electrodes should be stored so that they are isolated from the patient.
- For surgical procedures where the high-frequency current could flow through parts of the body having a relatively small cross-sectional area, the use of bipolar techniques may be desirable in order to avoid unwanted coagulation.
- The output power selected should be as low as possible for the intended purpose.
- The use of flammable anaesthetics or oxidizing gases such as nitrous oxide (N<sub>2</sub>O) and oxygen should be avoided if the surgical procedure is carried out in the thorax area or the head, unless these agents are sucked away.
- Non-flammable agents should be used for cleaning and disinfection wherever possible.
- Flammable agents used for cleaning or disinfecting, or adhesive solvents, should be allowed to evaporate before starting the electrosurgical procedure. There is a risk of pooling of flammable solutions under the patient, in body depressions (such as the umbilicus) or in body cavities (such as the vagina). Any fluid pooled in these areas should be dried up before using high-frequency electrosurgical equipment. Attention should be paid to the danger of ignition of endogenous gases. Some materials, for example cotton, wool and gauze, when saturated with oxygen may be ignited by sparks produced in normal use of the high-frequency electrosurgical equipment.
- For patients with cardiac pacemakers or other active implants, there is risk of damages or interferences with the pacemaker operation. In case of doubt, approved qualified advice should be obtained.
- Interference produced by the operation of high-frequency electrosurgical equipment may adversely influence the operation of other electronic equipment.
- Failure of the high frequency electrosurgical equipment could result in an unintended increase of output power.
- The qualified operator of the device should be a physician or medical workers under the supervision of the physician.
- Keep active accessories away from patient when not in use, never lay the pencil on the patient or the patient drape. Use a safety holster or place the pencil on the instrument tray when not in use.
- Keep active accessories away from flammable objects, gases or vapors at all times when in use.
- Keep desired voltage/power as low as possible to achieve the desired end effect in order to minimize the potential for capacitive coupling and inadvertent burning at high voltages.
- The cable be positioned in such a way that contact with the patient or other cable is avoided. Temporarily unused active electrodes should be stored in a location that is isolated from the patient.
- Accidental activation or use the equipment when the power supply exceeds the rated voltage, the patient, operator, or his assistant may get hurt, and the electrode may be damaged.
- Prevent soaking the product or the contact with big amounts of liquid, product functionality may be affected.

#### OPERATION ENVIRONMENT

Environmental temperature: 5 to 40 °C / 41°F to 104 °F  
 Relative humidity: ≤ 80% RH  
 Atmospheric pressure: 860 - 1060 hPa

#### STORAGE CONDITIONS

Temperature requirement in warehouse: -20°C to 60°C / -4°F to 120°F.  
 Humidity requirement in warehouse: 10 - 90 %RH.  
 Atmospheric pressure requirement in warehouse: 500 - 1060 hPa

#### DISPOSAL OF WASTE

The used disposable products should be destroyed and treated according to relevant local / national regulatory requirements. Randomly discard is strictly forbidden.

#### SHELF LIFE

The shelf-life of this disposable electrosurgical pencil can be guaranteed as 5 years within the established storage conditions.

#### EMC DECLARATION

This product should not be used adjacent to or stacked with other unit. If adjacent or stacked use is necessary, this product should be observed to verify normal operation in the configuration in which it will be used.

#### GUIDANCE AND MANUFACTURE'S DECLARATION – ELECTROMAGNETIC EMISSIONS

The Electrosurgical pencil is intended for use in conjunction with ESU suitable for use in the specified electromagnetic environment (s) and according to the following standard's emission requirements.

Phenomenon	Profession healthcare facility environment
Conducted and radiated RF emissions	Group 1, Class A
Harmonic distortion	IEC 61000-3-2, Class A
Voltage fluctuations and flicker	IEC 61000-3-3

#### GUIDANCE AND MANUFACTURE'S DECLARATION – ELECTROMAGNETIC IMMUNITY

The Electrosurgical pencil is intended for use in conjunction with ESU suitable for use in the specified electromagnetic environment (s) and according to the following immunity test levels. Higher immunity levels may cause the electrosurgical pencil and ESU's essential performance lost or degraded.

Phenomenon	Basic EMC standard or test method	Professional healthcare facility environment
Electrostatic discharge	IEC 61000-4-2	+/- 8 kV contact +/- 2 kV, +/- 4 kV, +/- 8 kV, +/- 15 kV air
Radiated RF EM fields	IEC 61000-4-3	3V/m 80MHz-2.7GHz 80%AM at 1kHz or 2Hz
Proximity fields from RF wireless communications equipment	IEC 61000-4-3	See the RF wireless communication equipment table in "Recommended minimum separation distances".
Rated power frequency magnetic fields	IEC 61000-4-8	
Electric fast transients bursts	IEC 61000-4-4	
		For input a.c. power port d.c. power lines or signal input/output lines whose length exceeding 3m
Surges	IEC 61000-4-5	
		For 1. input a.c. power port; 2. all d.c. power ports connected permanently to cables >3m 3. output signal output lines connected directly to outdoor cables
Conducted disturbances induced by RF fields	IEC 61000-4-6	
		For 1. input a.c. power port; 2. all d.c. power ports connected permanently to cables >3m 3. all patient-coupled cables 4. SIP/SOP whose maximum cable length ≥ 3m
Voltage dips	IEC 61000-4-11	
Voltage interruptions	IEC 61000-4-11	

UT: rated voltage(s); E.g. 25/30 cycles means 25 cycles at 50Hz or 30 cycles at 60Hz

#### RECOMMENDED MINIMUM SEPARATION DISTANCES

Nowadays, many RF wireless equipments have been used in various healthcare locations where medical equipment and/or systems are used. When they are used in close proximity to medical equipment and/or systems, the medical equipment and/or systems' basic safety and essential performance may be affected. This medical equipment and/or system has been tested with the immunity test level in the below table and meet the related requirements of IEC 60601-1-2:2014. The customer and/or user should help keep a minimum distance between RF wireless communications equipment and this medical equipment and/or systems as recommended below.

Test frequency (MHz)	Band (MHz)	Service	Modulation	Maximum power (W)	Distance (m)	Immunity test level (V/m)
385	380-390	TETRA 400	Pulse modulation 18Hz	1.8	0.3	27
450	430-470	GMRS 460, FRS 460	FM ± 5 kHz deviation 1 kHz sine	2	0.3	28
710	704-787	LTE Band 13, 17	Pulse modulation 217Hz	0.2	0.3	9
745						
780						
810	800-960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18Hz	2	0.3	28
870						
930						
1720	1700-1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation 217Hz	2	0.3	28
1845						
1970						
2450	2400-2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217Hz	2	0.3	28
5240	5100-5800	WLAN 802.11 a/n	Pulse modulation 217Hz	0.2	0.3	9
5500						

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TELIC, S.A.U.
Zhejiang Jinhua Huatong Medical Appliance Co., Ltd.
Shanghai International Holding Corp.
GmbH (Europe)