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Disposable Adhesive Surface Electrodes

Instructions For Use & Product Information

EN

Indications For Use

The Disposable Adhesive Surface Electrodes are intended for use with recording, monitoring and stimulation equipment for the purpose of stimulating/recording of biopotential signals. Electrodes are applied in the study of biopotentials such as Electroencephalograph (EEG), surface Electromyography (EMG), nerve conduction and Evoked Potential signals (EP). Electrodes are non-invasive as they are applied to the patients skin using a self-adhesive solid-gel surface. The electrodes are non-sterile and for single patient use only.

Caution/warnings

- In the USA Federal Law restricts this device to sale by or on the order of a physician or licensed practitioner.
- Only to be used by trained and skilled personnel, under responsibility of a physician or licensed practitioner.
- If irritation develops, discontinue use and consult your clinician.
- The electrodes are non sterile.
- The size, shape and type of electrode may affect the safety and effectiveness of electrical stimulation and recording.
- The electrical performance characteristics of electrodes may affect the safety and effectiveness of electrical stimulation and recording.
- Using stimulation electrodes that are too small or incorrectly applied could result in discomfort or skin burns.
- The user should contact the manufacturer of the electrical stimulator or recording device if the user does not know if the electrode can be used with the stimulation or recording device.
- Do not share electrodes with other persons because of the risk of adverse skin reactions and disease transmissions.
- Storage temperature between 15 °C(59 °F) and 35 °C(95 °F).

Instructions for application

- Apply only to unbroken skin. Thoroughly cleanse application site to remove lotion, dirt, etc. Dry site completely. Clip hair if necessary. **DO NOT SHAVE OR ABRASE AREA!**
- Remove electrodes from package. Remove electrodes from the liner by lifting the edge of the electrode.
- Apply electrodes to treatment site by firmly pressing down centre of pad first and then smoothing down to electrode edges.

Instructions for removal

- Turn off device.
- Disconnect lead wire from device.
- Remove electrode from skin.
- Remove possible residue of the gel on the skin with some water and a normal soap.

Technomed Europe
Amerikalaan 71
6199 AE Maastricht-Airport
The Netherlands
www.technomed.nl
T: +31 43 408 68 68