DORMO®-STRIP

	TELIC, S.A.U. Polígono Industrial Can Barri C/ Molí d'en Barri, 7 08415 Bigues i Riells Barcelona (Spain) www.telic.es
--	---

WITH THE GUARANTEE OF
TELIC GROUP

-MUB80/20221003-

ES. EXTRACTOR DE VENAS DE UN SOLO USO

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EN. SINGLE-USE VEIN STRIPPER

INDICATIONS

The vein strippers are a single-use medical device sterilized by ethylene oxide, intended to be used for surgical stripping of varicose veins.

This product is not intended to be used in central venous system.

Different models are available for use with the two most common stripping techniques: conventional stripping and invagination stripping.

REF	Description
VE-022	Conventional stripping
VE-022 OP	Conventional stripping
VE-025	Invagination stripping

INSTRUCTIONS FOR USE

Do not open the package until immediately prior to use.

From the set of components, take one of the cables, the handle and an adequate olive according to the size of the vein to be stripped.



The olive and the handle must be attached at opposite ends of the cable after the cable has been inserted in the vein to be removed.

Olive assembly: Insert the cable into the slot of the olive. A "click" must be heard to ensure that it has been correctly inserted. Seat the stopper, located at the end of the cable, into the olive so that it remains locked.

Handle assembly: Insert the cable into the slot of the handle. Seat the stopper, located at the end of the cable, into the handle.

For invagination stripping, instead of using the olive, the end of the vein to be removed must be knotted in the orifice of the stopper using a suture.

Pull the cable from the end where the handle is located to strip the vein.

Conventional stripping	
Invagination stripping	

TIME OF USE

Product for short-term use, intended for continuous use during the surgical intervention.

WARNINGS

This product must be used by authorised medical staff.

Follow the instructions for use and disposal of the product.

Do not use the product for surgical interventions of the central venous system.

Verify that the stopper located at the end of the cable is adjusted to the internal diameter of the olive.

Once stripping has started, do not push back on the cable.

There might be a hooking risk when introucing the cable into the vein. It is recommended to use another stripper to remove the first.

Product for single use only, sterilized by ethylene oxide. Do not re-use. Do not re-sterilize.

The re-use or re-sterilization of this product could cause:

- Cross infection due to the presence of biological waste;
- Alterations to materials;
- Loss of functionality of the product.

Re-sterilization of the product is not guaranteed.

Do not use if the individual packaging is open or damaged.

After having checked the integrity of the packaging, verify its content. If the packaging shows signs of damage or visible defects, do not use the product and return it to TELIC, S.A.U.

Check the expiry date on the packaging. Do not use after the expiry date.

The recommended storage conditions are between 5º C and 40º C (41º F – 104º F). The relative humidity restrictions are 20-80%.

DISPOSAL

– Dispose of the product following the hospital protocol for contaminated waste, in accordance with current standards.

– Dispose of the packaging following recyclable waste management protocols.

WARRANTY

TELIC, S.A.U. guarantees that the product complies with the Regulation (EU) 2017/745 on Medical devices and was manufactured following the procedures of the TELIC, S.A.U. quality management system certified under ISO standard 13485. TELIC, S.A.U. accepts no liability for any medical costs or direct or indirect harm due to the lack of operation or incorrect operation of the products when used without observing the instructions for use. We advise you to duly report any operational problem or product fault to the Quality Assurance Department of TELIC, S.A.U. In the event of any serious incident related to the product, this must be reported immediately to TELIC, S.A.U. (incidents@telic.es) and/or the competent authority of the Member State in which the incident occurs.

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