Member State in which the product was manufactured. This must be reported immediately to TELIC, S.A.U. (incident@telics.es) and/or the competent authority of the problem or product that to the Quality Assurance Department of TELIC, S.A.U. In the event of any serious incident related to the operation of the products when used without observing the instructions for use. We advise you to duly report any operational TELIC, S.A.U. accepts no liability for any medical costs or direct or indirect harm due to the lack of operation or incorrect manufactured following the procedures of the TELIC, S.A.U. quality management system certified under ISO standard 13485. TELIC, S.A.U. guarantees that the product complies with the Regulation (EU) 2017/745 on Medical devices and was GUARANTEE - NOTICE

- If they are stored in this way, the electrodes can be used safely up to 7 days after the bag has been opened.
- Once the bag has been opened, use the electrodes immediately.

STORAGE

ECG Electrodes may be used on patients for up to 24 hours.

DURATION OF USE

- During surgical procedures electrodes should be placed as far as possible from the electro- surgical area to minimize unwarmed Rf current flow and burnings in the patient.
- Electrodes may not be repositioned.
- In the event of use over an extended period, check for secure adhesion and absence of skin reaction on a daily basis.
- ECG monitoring devices. Never allow contact with an external electrical source.
- Use electrodes with no lead wires, only in combination with approved patient cables, adapters and
- Do not use the electrodes if the pouch is found damaged.

PRECAUTIONS

- Do not use if the gel is dry or discoloured.
- Dissolve the adhesive.
- Careless removal of electrode may cause damage to the skin. Should the electrode be difficult to remove, use water to Gently removal of electrodes from the pouch to be ensured so as to prevent any injury to the fingers as plastic cover is fine cut.
- Compatible with conventional cleaning and sterilization procedures.
- Procedures may leave harmful residues or cause malfunction of the device. The design and materials used are not For single use only. Use on other patients can cause cross infection. Do not soak, rinse, or sterilize this device, as these Plastic cover should be immediately disposed off in the assigned waste bin for disposal after removal of electrodes.
- The product is not for internal use. In case of accidental ingestion of the electrode consult a physician immediately.
- For use only on intact, clean skin (do not use on open wounds, lesions, infected or inflamed areas).
- Any allergic or other skin reactions should be immediately evaluated by a physician.
- ECG electrodes should only be used by or in consultation with a health care provider familiar with their proper placement

WARNINGS

- Recycle package.
- Dispose off as hospital/ contaminated waste, in accordance with regulations in force.

DISPOSAL

- Remove gel residue, if required, with a paper towel or water.
- After monitoring, lift the tab to gently remove the electrode from the skin. ⑤
- Turn off equipment on and perform the ECG procedure. ④
- Place the electrodes on the skin, smooth the backing surface with the fingers as necessary. ③
- Connect the lead wires to the electrodes and peel the plastic liner. ②
- application of the electrodes. ①
- If necessary, clean the skin with appropriate preparation material. Always ensure the application site is dry prior to fluid, dirt or other impediments to good electric contact.
- Determine the application sites. Remove any excessive hair, calluses and ensure that the skin is clear of any oils, creams, lotions, dirt or other impediments to good electric contact.
- Open the pouch along the indicated line.

APPLICATION

- suitable for use on newborn or premature infants.
- Due to high adhesion strength of the adhesive material, this product is not suitable for use on patients with fragile skin. Not

AREA FOR APPLICATION / CONTRAINDICATIONS

- For use in passive sensing of bio-electric signals from the intact skin.

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