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## EN. CONDUCTIVE GEL FOR ELECTRODES

Conductive gel for electro-medical applications (ECG, TENS).

The product is marketed with the format and reference below.

| REF    | Description                       |
|--------|-----------------------------------|
| 602047 | 250 ml bottle with dispensing cap |

### INSTRUCTIONS FOR USE

1. Clean and dry the area of application, if necessary.
2. Apply a small amount of gel to the conductive area of the electrode.
3. Connect the electrodes to the device to carry out the procedure as per the usual protocols.
4. Remove the electrodes, wipe off the gel with a single-use cloth then dispose of the product in accordance with the corresponding protocol. If the electrodes can be reused, follow the cleaning instructions given by the manufacturer.

### TIME FOR USE

The electrodes should be replaced if a decrease in quality is detected or at the latest after 24 hours.

### WARNINGS

The gel should only be used by or under the supervision of a healthcare professional who knows how to apply and properly use it. – Do not use on pre-gelled electrodes. – If an allergy or any other skin reaction occurs, it must be immediately evaluated by a doctor. – Follow the instructions for use and the instructions for disposal of the product. – For use on healthy and clean skin only (do not use on wounds, injuries, infected areas or skin inflammations). – The product is not intended for internal use. – Consult a doctor immediately in case of accidental ingestion. – Single-use gel. Reuse on different patients can result in cross infections. – Do not immerse, rinse or sterilise this product, as these procedures can leave harmful residues or cause a dysfunction. Neither the design of the product nor the materials with which it is manufactured are compatible with conventional cleaning and sterilisation procedures. – Do not use the product if the gel is dry or has changed colour. – The recommended temperature limits for its storage are between 5°C (41°F) and 40°C (104°F) – Check the expiry date printed on the packaging. – Improper use or incorrect application of the gel can result in incidents and/or make the treatment ineffective.

### DISPOSAL

The product must be disposed of in accordance with current standards. Dispose of the packaging as per the regulations on the management of recyclable waste.

### WARRANTY

TELIC, S.A.U. guarantees that the product complies with the Regulation (EU) 2017/745 on Medical devices and was manufactured following the procedures of the TELIC, S.A.U. quality management system certified under ISO standard 13485. TELIC, S.A.U. accepts no liability for any medical costs or direct or indirect harm due to the lack of operation or incorrect operation of the products when used without observing the instructions for use. We advise you to duly report any operational problem or product fault to the Quality Assurance Department of TELIC, S.A.U. (incidents@telic.es). In the event of any serious incident related to the product, this must be reported immediately to TELIC, S.A.U. and/or the competent authority of the Member State in which the incident occurs.

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