# blayco®

ES. PLACA ELECTROQUIRÚRGICA PREGELADA

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## EN. PRE-GELLED ELECTROSURGICAL PLATE

### INDICATIONS

Electrosurgical plates are non-sterile medical devices for single use only whose intended use is to act as closing element in the electrosurgical circuit. The plates collect the electrosurgical current delivered to the patient and returns it to the electrosurgical generator.

The plate provides a large contact surface with the patient, compared with the active electrode, that allows reducing the current flow density and minimize the risk of electrosurgical effects or burnings.

The product is commercialized under different models that depend on the patient or electrosurgical procedure. Dual plates are compatible with the contact quality monitoring system (REM or ARM).

All plate models are available with tab connection or with a single use cable. For references with cable, several connector models are available, compatible with the main electrosurgical units available in the market.

### INSTRUCTIONS FOR USE


Remove the plate from the package just immediately before use. Select a well vascularised area near the incision. Do not place on bony protrusions, metallic prosthesis, tattoos or scars. Avoid if possible the areas with a lot of fat under the skin. Shave if necessary. Clean and dry the application area. In the case of the elderly, smooth out the skin beforehand.

Check the adequate conditions of the plate, the cable and the plug. Remove the protective liner and check that the hydrogel is in good conditions before use. Apply the adhesive face of the plate smoothing and pressing lightly.

Connect the plate to the electrosurgical generator using the corresponding cable. Check that electrode has been conveniently detected by the electrosurgical equipment.

Proceed with the electrosurgical intervention according to usual medical protocols. Once the intervention is finished, remove the plate gently.

### ENVELOPE OPENING

1. Tear the envelope through vertical welding
2. Tear the upper horizontal welding and open the envelope
3. Extract the plate

### PLATE APPLICATION

1. Shave the skin area where the plate will be applied, if its necessary
2. Remove protective film from the plate
3. Apply plate on selected area
4. Connect to electrosurgical unit

### DURATION OF USE

Product for short-term use, intended for continuous use during the electrosurgical intervention.

### WARNINGS

This product must be connected and used by qualified medical staff. For using the electrosurgical generator, pencils, cables and other accessories, please refer to the instructions for use supplied by the manufacturer. Follow the instructions for use and disposal of the product. Do not open the package until immediately prior to use to avoid gel drying.

Check the adequate conditions of the plate, the cable and the plug. Avoid placing the plate in areas of skin over metallic implants, scar tissue, adipose tissue, bony prominences, tattoos, erythema and lesions of any kind.

Do not apply the plate in areas where fluids may access. Do not cut or modify the plate in any case. Prevent the active cable from coming into contact with the patient's skin or becoming entangled. Do not activate the electrosurgical equipment until the plate is well adhered to the patient. Do not use electrode gel.

If the patient is repositioned, verify the contact of the plate and the cable connections. Place the E.C.G. electrodes as far as possible and equidistant from the incision area. During the surgical procedure always try to select the lowest energy level that will produce the desired effect. If power adjustments greater than usual are required, first check the plate-patient contact and inspect clamps, cables, connections and active accessories. Do not increase the high-frequency output before making the following checks:

- Proper positioning of the neutral plate;
  - Proper insertion of the cables and their connectors;
  - Proper activation of the operating keys (manual or foot control);
  - That there is no damage to the cable insulation;
  - That the electrode and the extension (if any) are not dirty.
- Always select the most adequate plate model for each procedure or patient.

Do not use electrosurgical plates in Magnetic Resonance Imaging facilities. Remove the plate before introducing the patient into Magnetic Resonance Imaging facilities.

Product for single use only. Do not reuse. Non-sterile product. The sterilization of this device is not guaranteed.

The reuse or sterilization of this device could result in:

- Possible presence of biological residue that could cause cross-

infections;

- Alteration of the materials;

- Loss of the product's functional features.

Do not use if the individual package is open or damaged.

After checking the integrity of the package, please check its content. In case of visible damage or defects, do not use the product and return it to TELIC, S.A.U.

Check the expiry date on the package. Do not use after the expiry date. The recommended conditions of use and storage temperatures are between 5°C and 35°C (41°F - 95°F). The relative humidity restrictions are 20-80%.

Pacemakers may be damaged by the electrosurgical current. Do not expose a patient with a pacemaker to an electrosurgical current without first consulting a cardiologist.

## ELECTROMAGNETIC COMPATIBILITY

Electrosurgical plates are intended for use in operating theatres of healthcare facilities equipped with high-frequency electrosurgical units. The product is manufactured upon some metallic components, then the use of the product is not compatible with Magnetic Resonance Imaging.

Intended use of the product is to safely remove electrosurgical current from the patient and return it to the generator. Performance of electrosurgical plates is not expected to be affected by electromagnetic disturbances. However, interferences in the electrosurgical system cannot be discarded and needs to be taken into account in case of observing failures in plates detection, unusual current settings or excessive heating of the application site.

Use of this electrosurgical equipment adjacent or stacked with other equipment or cables of other equipment should be avoided because it could result in improper operation because of electromagnetic disturbances. If such use is necessary, the equipment should be observed to verify that they are operating normally.

The longest cable supported by electromagnetic disturbance testing is of 5m. For references with tab connection, intended to be connected to the electrosurgical generator with a reusable clamp cable, do not use cables longer than 5m. Use of accessories or cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions and result in improper operation.

Portable radio-frequency communications equipment (including peripherals such as antenna cables and external antennas) should be used not closer than 30 cm (12 inches) to any part of the electrosurgical system, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

The emissions characteristics of this equipment make it suitable for use in industrial areas and hospitals (EN 55011/CISPR 11 class A). If it is used in a residential environment (for which EN 55011 / CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such

as relocating or re-orienting the equipment.

The product is classified as group 1, class A product according to EN 55011/CISPR 11. This means that this product does not generate and/or use intentionally radio-frequency energy, in the form of electromagnetic radiation, inductive and/or capacitive coupling, for any purpose and that it is suitable for use in all establishments other than domestic and those directly connected to a low voltage power supply network which supplies buildings used for domestic purposes.

There may be potential difficulties in ensuring electromagnetic compatibility in environments other than hospital or industrial, due to conducted as well as radiated disturbances.

For product construction materials and performance, it is not expected that electrosurgical plates are degraded due to electromagnetic disturbances. Accordingly, electromagnetic immunity tests are not required for this product.

Operation of electrosurgical plates can be guaranteed into the following range of frequency output: 200 kHz – 5.000 kHz.

## DISPOSAL

- The product must be disposed of according to hospital protocols for contaminated waste, in accordance with regulations in force.

- Dispose of the packaging according to recyclable waste management.

## WARRANTY

TELIC, S.A.U. guarantees that the product complies with the Regulation (EU) 2017/745 on Medical devices and was manufactured following the procedures of the TELIC, S.A.U. quality management system certified under ISO standard 13485. TELIC, S.A.U. accepts no liability for any medical costs or direct or indirect harm due to the lack of operation or incorrect operation of the products when used without observing the instructions for use. We advise you to duly report any operational problem or product fault to the Quality Assurance Department of TELIC, S.A.U. In the event of any serious incident related to the product, this must be reported immediately to TELIC, S.A.U. (incidents@telic.es) and/or the competent authority of the Member State in which the incident occurs.

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