REF NEK300A LOT 1034/12-22 2027-12-11

(01) 8427734038109 (17) 271211 (10) 1034/12-22

WITH THE GUARANTEE OF TELIC GROUP

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### EN. REUSABLE CABLES FOR ELECTROSURGICAL PLATES

TELIC, S.A.U. markets different models of reusable cables, which include different connectors compatible with the main electrosurgical units available on the market.

All references include a clamp for connection to the tab of the return electrode. The cables are available in lengths of 3m or 5m according to the following references:

Connector	Equipment	REF 3m	REF 5m
	ERBE type connector.	NEK300D	NEK500D
	VALLEYLAB type connector. Dual plates.	NEK300A	NEK500A

#### INSTRUCTIONS FOR USE

- 1. Inspect the clamp, cable and connector before using. Check that there are no cuts or corrosion residue.
- 2. Apply the plate to the patient in accordance with the instructions.
- 3. Position the lever in the release position and insert the tab of the plate between the jaws of the clamp, so it is well centred.
- 4. Place the lever in the desired position.
- 5. Verify the alarm circuit of the electrosurgical unit works correctly by starting up the unit before connecting the cable to it.
- 6. Check that the alarm stops when the cable is connected. If the alarm does not stop, check the plate-patient contact, cable, connections and generator.

#### TIME FOR USE

As with the electrosurgical plates, it is not advisable to use the product continuously for more than 24 hours.

#### CLEANING PROCEDURE

- 1. Clean with a brush or sponge dampened with alcohol solution or soap and water. Do not use aggressive chemical products.
- 2. Rinse off the soap and water with distilled or deionised water.
- 3. Ensure that the clamp and the cable are perfectly dry after cleaning. During the cleaning procedure, do not allow the liquids to come into contact with the connector.

#### WARNINGS

The reusable cables must only be used by a health professional familiar with its proper use and application, or under their direction. - Verify the connector type and that the cable is in good condition before use. Do not use if the product is found to be damaged. - Do not cut or modify the cable under any circumstances. - Do not allow the cable to come into direct contact with the patient's skin or become tangled. - Do not operate the electrosurgical unit until the plate is properly adhered to the patient. - If the patient is repositioned, check the plate contact and the cable connections. - In patients with pacemakers, consult the cardiologist before subjecting the patient to an electrosurgical current. - During the surgical procedure, always try to select the lowest energy level that gives the desired effect. - If power has to be adjusted above the usual level, first check the plate-patient contact and inspect the clamp, cables, connections and accessories in use. - Once used, dispose of in accordance with the current regulations on medical devices. - Store in a cool, dry place. - Check the expiry date on the container.

#### DISPOSAL

- Dispose of the product if any indications of damage or malfunction are observed.
- Dispose of the product following the hospital protocol for contaminated waste, in accordance with the current standards.
- Dispose of the packaging following the recyclable waste management protocols.

#### WARRANTY

TELIC, S.A.U. guarantees that the product complies with the Regulation (EU) 2017/745 on Medical devices and was manufactured following the procedures of the TELIC, S.A.U. quality management system certified under ISO standard 13485. TELIC, S.A.U. accepts no liability for any medical costs or direct or indirect harm due to the lack of operation or incorrect operation of the products when used without observing the instructions for use. We advise you to duly report any operational problem or product fault to the Quality Assurance Department of TELIC, S.A.U. (incidents@telic.es). In the event of any serious incident related to the product, this must be reported immediately to TELIC, S.A.U. and/or the competent authority of the Member State in which the incident occurs.

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