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EN. DISPOSABLE ELECTRODE TIP CLEANER
INDICATIONS Disposable abrasive sterile pad intended for electrosurgical pencil cleaning during an electrosurgery intervention. The product is sold under the following models and references:

REF	Description
AL-40	50 x 50 mm abrasive pad

INSTRUCTIONS FOR USE

- Do not open the packaging until immediately prior to use.
- Remove the protective film from the adhesive back side of the pad and press the tip cleaner onto the operation drapes so that it is firmly in place to facilitate cleaning of the electrosurgical pencil tip.
- NB: To prevent eliminated material from falling into the operating field, do not perform tip cleaning in the immediate vicinity of the operating area.
- The radiopaque tracer inside the product allows it to be located by X-ray in case of accidental loss within the operating field.
- If a reduction in the abrasive and cleaning potential of the pad is observed during use, it is advisable to replace the product.
- The product must be disposed after the surgical intervention.

DURATION OF USE

Product for short-term use, intended for continuous use during the electrosurgical intervention.

WARNINGS

This product must be used by qualified medical staff.

For using the electrosurgical generator, pencils, neutral plates, cables and other accessories, please refer to the instructions for use supplied by the manufacturer.

Follow the instructions for use and disposal of the product.

Product for single use only, sterilized by ethylene oxide. Do not reuse. Do not re-sterilize.

The reuse or re-sterilization of this device could result in:

- Possible presence of biological residue that could cause cross-infections;
- Alteration of the materials;
- Loss of the product's functional features.

Re-sterilization of the product is not guaranteed.

Do not use if the individual package is open or damaged.

After checking the integrity of the package, please check its content. In case of visible damage or defects, do not use the product and return it to TELIC, S.A.U.

Check the expiry date on the package. Do not use after the expiry date.

The recommended conditions of use and storage temperatures are between 5° C and 40° C (41° F - 104° F). The relative humidity restrictions are 20-80%.

Keep the active electrode tip in light contact with tissue. Do not apply mechanical force to the electrode.

DISPOSAL

- The product must be disposed of according to hospital protocols for contaminated waste, in accordance with regulations in force.
- Dispose of the packaging according to recyclable waste management.

WARRANTY

TELIC, S.A.U. guarantees that the product complies with the Regulation (EU) 2017/745 on Medical devices and was manufactured following the procedures of the TELIC, S.A.U. quality management system certified under ISO standard 13485. TELIC, S.A.U. accepts no liability for any medical costs or direct or indirect harm due to the lack of operation or incorrect operation of the products when used without observing the instructions for use. We advise you to duly report any operational problem or product fault to the Quality Assurance Department of TELIC, S.A.U. In the event of any serious incident related to the product, this must be reported immediately to TELIC, S.A.U. (incidents@telic.es) and/or the competent authority of the Member State in which the incident occurs.

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